

Product datasheet for **SC205286**

C9orf174 (CCDC180) (NM_020893) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C9orf174 (CCDC180) (NM_020893) Human 3' UTR Clone
Symbol:	C9orf174
Synonyms:	C9orf174; CFAP76; FAP76
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_020893
Insert Size:	1887 bp



[View online »](#)

Insert Sequence: >SC205286 3'UTR clone of NM_020893
The sequence shown below is from the reference sequence of NM_020893. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTGCACACTATCCAAGGCTGTATGTGTACCCTCCGCCCCACCATGAATAAACACTTTCTTATACAGA
CTCCTTCCCTGTCCATCTACCTGCCTACCTACTTTCCGTCCATCCCTCCCTTCATCCCTCCACCC
CTGCTCTGTGCCGGGCACTGTAGCTTTACCAGCGAACAGGACACAGCATGGTCCCTGCCACGTGGAGC
CCTCTTCCCATGAGGAAGGCAAGCATGAAAGGGGAATTCCACGTGGTTAGCAGGGAGGGTGAAGTACCAG
GAATCAGTCACCTTCAAATAGGCACAGAGCTTCAAAGGACATCACCTCCAACAGCATGCTGGGCAGGG
CACGTCATCTTGAAAAATGAAATAAAGGTCATGGGAAAAATGATATTTTAATTTAACAGAATAACAAGT
ATGGGTTATAATTGCAACTCTGCACTGGTTCCTAGACCCAACAAATAACTAATCCAATGTTGATTGGTA
TTTGTGCCACAATTTGTCTGAGCAAATTGCAATATTTTGGAAAAAATAATCTCCAGGCTCCTTAG
TCTGTGTGTCAGATCATTCTCTGGAGTCTCCAGCCTGCCAGTCTCTTGTTCGTCAACTCTCTTC
TTCAATTACATGGCGTGAGCGTACAGCCTTCATTTGTCGGTAGCTGTGTGTGTCTCTGGACAGTTTC
CTCAGCCATAGACCTAGGAGGCCCGAGCTCCCAATTTCCAGCCCTGAAAGCTGTGTGACCCTCCAAG
TGCAGTTCCTCACCTGGGGCACTCACTCTGCTCTCCTGGTTTGTCCACCCAGTGAGGCCCTGCTTG
GTTAGCTTGGTTAGCTTGGTGGCCTTTGTGCACTGAAGCCAGGAGAGTGAGGAGTCCCCTTCCATCAAG
GCAAAAGCAATCTCCTGTTCTTTCTGTGCTCCCCCACCCCAACCAATGTGTACCTTTCATTGCT
TATTATCCGAAAGGAATATATGCGCATTGTAGAAAAAGAATTTATGAAACCAGCAAAAGGATAACTT
TAAATCCTGAGCCGGGCGGGTGGCTCACGCCTGTAATCCTAGCACTTTGGGAGGCCAAGGCGAGCAG
ATCATCTGAGCTCAGGAGTTCGGGGCCAGCTGAGGCAACATGACGAAACCCATATCTACTAAAAATAC
AAAAAGTTAGCCAGGCATGGTGGCACACCTGTAATCCAGCTACTGGGAGGCTGAGGCATGAGAAT
CGCTTTAACCTGGGAGGTGGAGTTGCACTGAGCCAAGATCATGCCACTTCACTCTGGCCTGAGCAACA
GAATGAGACTCTGTCTCAAATAAACGAATAAATAAATAAATAAATCCTGGAATCTAGAGACAGTCACTG
TCCATAGCATCTTCCATGTATGTGAGTTATTTTATAAAATGGAATCATAACATGCATTAGGTTCCAT
AACTGCTTGGTAAATGTTAACTGGATAAACCTATGCCATATCCTTCAGTATTCATTTACATCATTGG
TTTTCACTTTTCTGACTCTGACCCACTGTATAAACATGTTCTAAACACAGGCTCCTTGACACATAC
ACAACCTGAGACCAACAGGTACAGGACAGCACTTAACCTTACTAGAAAGGCTGCTTTCTAATATTTTC
TCCTGCTTTTTAAATCCTAGCCAGTCCACTAACTCAGTTCATGGCCCAAAAGGGCTTTGATTGCA
GTGTGGGACACTGTCAATTAAACAGCTGCAGAGTGTCCCACTGCCTGGGGAACAGATCGCCTCTGTC
GGACAAGGGCGAGGATGTCCATGCAGGCTCCTAGAAAGAAGGCAGAGATGACATCCTCCACATTAG
CTCTTTGTACACATCTGTAATCAT
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_020893.6](#)

Summary: The protein encoded by this gene contains a coiled-coil domain. Alternative splicing results in multiple transcript variants encoding different isoforms. A single nucleotide polymorphism (SNP) in this gene has been associated with increased susceptibility to Behcet's Disease (PMID: 19442274). [provided by RefSeq, Dec 2016]

Locus ID: 100499483

MW: 69.9