

Product datasheet for SC205290

ORF1 FL49 (CYSTM1) (NM_032412) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ORF1 FL49 (CYSTM1) (NM_032412) Human 3' UTR Clone
Symbol:	ORF1 FL49
Synonyms:	C5orf32; ORF1-FL49
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_032412
Insert Size:	408 bp
Insert Sequence:	>SC205290 3'UTR clone of NM_032412 The sequence shown below is from the reference sequence of NM_032412. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TGCTGCTGTCTCTGGGACATGCTCACCTGACCAGACCAGCCAGCCGTCCTGCTCCTGCCAGCTCTGCTG
CCACCTCTGACAGGTGTGCCTGCCCCATCTCTTCTGATTGCTGTTAACAAATGACTAGCTTTGCACAG
ACACCTCTACCTTCAGCACTATGGGATTCTAGATTAATGGGGTTGCTACTGTTTAATTCAGTGACTTG
ATCTTTTTTAATGTCCAAATCCATTTCTTATTGATCTTTAAAGATGTGCTAAATGACTTTTTTGGCCAA
AGGCTTAGTTGTGAAAAATATAATTTTTAAATTATACATTCAAGGTAGTGGCCAAATGTAACACATCAA
TCATGGAATGATTCTCTGCTAACAGCCGCTGTATGTTTCAATAAATTTGTCCAAAGCTCAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

RefSeq: NM_032412.4
Locus ID: 84418
MW: 15.4