

Product datasheet for **SC205297**

SPSB3 (NM_080861) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: SPSB3

Synonyms: C16orf31; SSB3

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_080861

Insert Size: 415 bp

Insert Sequence: >SC205297 3'UTR clone of NM_080861
The sequence shown below is from the reference sequence of NM_080861. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCCGAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TGCCAGAGGAAGCGCTGCCGCCGACCTGACTGACTTCCCAGTGGAAGTGCCTTCTTGGGCTGGGACAG
CCCCTTTCTCTGTCCCTTCTTTCTGTCCCTTCCTCCAGCCACACTCCAGGGCGGAGTTGGATGAG
GCCCGTCCGGAGGGAGCCATCTCTTGTCCCGAGGCTGGGACAGTCCTTTCTGTGGGGGCTCTAGGGCC
CCTCTGCTGTGTGCTGGGTGGGAAGCGGCTGCCCTGAGCCCCAGGTCTTGTGGGAGGCTGCGAGGAC
GAGAGCCTGGCTGGAGCCGCGTTGTGTTCCACAGGGCCTCGGTTTTTCTAACTTGCTCTGCATGC
TGTCAGCGGCTGCCCCGCCGTCATAGACTTAAAGGACTGCAATAAATGTAGAGTTGATGTCTAACACCCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_080861.4</u>
Summary:	May be a substrate recognition component of a SCF-like ECS (Elongin BC-CUL2/5-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.[UniProtKB/Swiss-Prot Function]
Locus ID:	90864
MW:	15