

Product datasheet for **SC205304**

LILRB5 (NM_001081443) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	LILRB5 (NM_001081443) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	LILRB5
Synonyms:	CD85C; LIR-8; LIR8
ACCN:	NM_001081443
Insert Size:	1395 bp



[View online »](#)

- Insert Sequence:** >SC205304 3'UTR clone of NM_001081443
The sequence shown below is from the reference sequence of NM_001081443. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site
- ```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGCATCTACGCCCCCTGGCCATCCACAGCCACGGGGACCCAGATCTCATACTCAACAGAAGGAGA
CTCAGAGACTCCAGAAGGCACAGGAGCTGCCCCAGTGGACACCAATGAACCCAGCCAGCCTGGACCC
CTAACAAAGACCACCAGGACATCCTGGAACTCTGGGACTCACTAGATTCTGCAGTCAAAGATGACTAA
TATCCTTGCATTTTTGAAATGAAGCCACAGACTTCTCAATAAATCAATGAGCTGAGAAAAGTAAACAG
AAATTAGAGCATGGTATAAATTTGGAATGATAATGTAATATTACACATTAATGATGAAATCGGAAAA
CTACAAATGAGCGAATGAATTAGAAAAGAATAAAACCTACGTAATTAATGACCTTGGCAATGACAGAAA
GAATTTAGAAAAAGAACAACAAATTATTCCAAATGAAGGTGTGAGGAAGGGGACAAAAATAACAAGAGG
AGTTACTAATGAGGGCTACGTGAAAACCTGATGAAGCCAAAAAGCTCATTCTTGAGAATGTGAATTAC
ATTCACAAATCCTAGCCACAATAAGCAAGGAAAAAAGCGGGTTCAGGCACACATTTCCATATGGGGGT
GAAACAGCAGACACCACCACAAATCTGACACATATTGCCTTTATTTTTTCACTTTTAAGTTCAGGGAT
ACATGTGCAGGTTTGTAGACAGATAAAGTGTCAAGGGGATTTGTCTTGGTTTTTGTGTGAGGGTT
TTTGTGTTTTGTTTTGTTTTGTTTTGAGACGGAGTCTCGCTCTATCACCAGGCTGGAGTACAGTGGAGT
GATCTCGGCTCACTGCAACCTCTGCCTCCCGGTTCAAGCGATTCTCCTGCCTCAGTCTCCCGAGTAAC
TGGGACTACAGGCACCTGCCACCACGCCCGCTAATTTTTGTATTTTAGTAGAGACAGGGTTTCTCCA
GTTTGGTCAGGCTGGTCTCAAACCTCCGACCTCAGGTGATCCGCCCGCTCAGCTCCCAAAGTGTGG
GCTTACAGGCGTGAGCCACCACGTCCAGCCATACATTTCAATTTTAAAGGGATGCGCCCTAGTCCTTA
GTTAGTCTCTCCTCATCTATAAAATGTTAGCTACTACCTCTTGGGCTATTGCTAGACATCGTTTTT
CTCTTCTTTCTTCTGACGCTACAATAGATAGGACATTCCTCCTCATTCTATTCTCCCAAGTACT
TTAAATTGCAATTTATAAAGTTTCTATGTACTCTAAAAAAATTCTGTTTTGTTTTCTAATTTTCAAT
AATTGGTGCTTCACTGTGCTTGTCTCGAAGGAATGAGTATTTGATTGTTCATTAATCTGATTT
TTCTATGTCTTCTAA
ACGCGTAAAGCGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```
- Restriction Sites:** SgfI-MluI
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
- Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
- RefSeq:** [NM\\_001081443.3](#)
- Summary:** This gene is a member of the leukocyte immunoglobulin-like receptor (LIR) family, which is found in a gene cluster at chromosomal region 19q13.4. The encoded protein belongs to the subfamily B class of LIR receptors which contain two or four extracellular immunoglobulin domains, a transmembrane domain, and two to four cytoplasmic immunoreceptor tyrosine-based inhibitory motifs (ITIMs). Several other LIR subfamily B receptors are expressed on immune cells where they bind to MHC class I molecules on antigen-presenting cells and inhibit stimulation of an immune response. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Locus ID: 10990  
MW: 54.3