

Product datasheet for SC205314

PRR24 (INAFM1) (NM_178511) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PRR24 (INAFM1) (NM_178511) Human 3' UTR Clone
Symbol:	PRR24
Synonyms:	PRR24
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_178511
Insert Size:	405 bp
Insert Sequence:	>SC205314 3'UTR clone of NM_178511 The sequence shown below is from the reference sequence of NM_178511. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GAGGCCGCGGAGGGGCGAAGACCCGGGTAACTCTCCCTTCCACCCCAACCCGGATCGCCAGCCCTCGAG
AGCTCTGTGCTCCACGCCGAGGATGCACCGTCTCTGGATTGGTCCGGCCTTCTTCTAATGACATCGCT
CAGCGTCTCTGGAGCCGTCATCCCGGGAATGGGGGCCAGGGGGTGTGAGCTCGGGGCTTGCCTCTT
GCAGCTACTCTGTGGTCAGGCCGGGTCTCCACCATCAGGAAGATCCCATCCTGAGCTCTGTCTCCTGC
CCCTCCTGCTGTGGGATGCTGAGCACAGAGCCACAGCCCATCTGCCTCTTACCTCCCTGAATCCGTG
TCCATCTGCAATAAACGACAGCCTCGGCTGCCTCGTGCTGTGTCTGTCTGCTGTTGTCTGCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

RefSeq: NM_178511.6
Locus ID: 255783
MW: 14.4