

Product datasheet for **SC205343**

TRIM73 (NM_198924) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: TRIM73

Synonyms: TRIM50B

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_198924

Insert Size: 427 bp

Insert Sequence: >SC205343 3'UTR clone of NM_198924
The sequence shown below is from the reference sequence of NM_198924. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TGGAAGTTCCACTCCATGGCCTCCAGGTAATTAACCTTGGAGAGAGCTCAGCCAGGGTCTGGTGGCTGCG
GGCAGGGGCATCTCAGCTCCACTGGTTCCTCCATTAGCTTAACCAGCGCTCCCAAGCAGCTGCCTAT
AGCTGGCTCTATACTGAGCTGGGGAAGATAGAGGAAAGTCACGTCCCTGCCTCAAGGGTCTCGCAG
ACAGGTGGGGAGGCAGATGGTGAAGTGTGGGTACCTAGAACAGCAGAAAGTCACTCAAGCTACAGAAAT
ACTAGAGGAGGGTAGCTCATGCCTGCAATCCCAGTACTTTGGGAGGCCAAGGCAGGAGTATTGCTGGAG
GCCGGGAGTTCGAGACCAGCCTGGCCAATGTAGTAACACCCCGTCTCTACAAAAATACAAAAATAAA
AAAATTAGTTGGG
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_198924.4</u>
Locus ID:	375593
MW:	15.8