

Product datasheet for **SC205357**

MINDY4 (NM_032222) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	MINDY4
Synonyms:	AQP-1; AQP1; C7orf67; CHIP28; FAM188B
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_032222
Insert Size:	412 bp
Insert Sequence:	>SC205357 3'UTR clone of NM_032222 The sequence shown below is from the reference sequence of NM_032222. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAACGCGATCGCC AACTGGAACGGCTCAGACCCCATCCTGTGACCGTTGGATGTGGGTAAACCCTGTGGTCCACCACTCATC ACCTCATCACCGAGGATGACAGCTGAACCCCAAGCCTCTGGGGCAGGTCTCATGTACCCCAACCTGGGT CAGCATGACTGCAGAAGCATCCAGAGCCTCCCTGCCCTTCCATGAAGGGCCCAACCAAGACTGTGCTG GGGACCCAGTGTGTTGCTGGGTCCCTCCAGCTGAGCTGTGACTGCTGAGTACTGGAAGGAGGTTGCC AGGGTCTCTGCTACCTTTGTCTGCATCCCTCCCTTGCTCCCTGCTGGGTGGTCCCTACCCAGGCCTCC AATGTGGTTGGCCTGGGGACTCATTACTTCTGGTAGCTGGCTTTCTATAAAATGGGTCAAACCA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_032222.3</u>
Summary:	Probable hydrolase that can remove 'Lys-48'-linked conjugated ubiquitin from proteins. [UniProtKB/Swiss-Prot Function]
Locus ID:	84182
MW:	14.3