

Product datasheet for **SC205420**

SLIT3 (NM_003062) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SLIT3 (NM_003062) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	SLIT3
Synonyms:	MEGF5; SLIL2; Slit-3; SLIT1; slit2
ACCN:	NM_003062
Insert Size:	2000 bp



[View online »](#)

Summary: The protein encoded by this gene is secreted, likely interacting with roundabout homolog receptors to effect cell migration. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2012]

Locus ID: 6586

MW: 75