

Product datasheet for **SC205435**

ZNF85 (NM_003429) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: ZNF85

Synonyms: HPF4; HTF1

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_003429

Insert Size: 425 bp

Insert Sequence: >SC205435 3'UTR clone of NM_003429
The sequence shown below is from the reference sequence of NM_003429. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
ATTCATACCGGAGAAAAATTACAAATATGAAAATTATGGCAAAGCTTTAATCAATTTACAAGTCTTACT
AAACATAAGAAAAATTTACTGGAGAGAACTACTAACCTGAAAGATGTGACAATAATTTTGACAACAC
CTCAGACTTATAAAAGTAATCATACTGGTGAGAAATCTAAAAATGTGAAGACTATGGCAAAGTCTTTA
AATGGTTGTCACACTTTAGGTAAGATAATTCATATTGGAACAACTACAAGTGCAACAATGTGGCAAA
ACTTAATTTATGCTCACACCTTACTGCACAGAAAAGATTTTGTAGTTGAGAAAAAGTATACAAATATAA
AGAATGTGGAAAAGCCGTTAATATCTGCTCACATCTTACTCAGCATCAGAAAGTACTTAATAAAGCAT
TATAAATGGAA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_003429.5</u>
Summary:	May be a transcriptional repressor.[UniProtKB/Swiss-Prot Function]
Locus ID:	7639
MW:	16.7