

## Product datasheet for **SC205438**

### **TMSBI5A (NM\_021992) Human 3' UTR Clone**

#### Product data:

Product Type:	3' UTR Clones
Symbol:	TMSBI5A
Synonyms:	Tb15; TbNB; TMSBI5; TMSBI5B; TMSL8; TMSNB
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_021992
Insert Size:	435 bp
Insert Sequence:	<p>&gt;SC205438 3'UTR clone of NM_021992</p> <p>The sequence shown below is from the reference sequence of NM_021992. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC CAAGAGAAAGAGTGTGTTCAAACATCAATAATGGGGATCGCCTCCCAACAGCAGATTTTCGACATTACC TGAGAGTCTTGATTTTAGGCTTGTTTTGTAAACCCATGTGTTGTAGAGATTTTAGGCGTCTTCGGA TATCTTCTCACCTATGTTCCCTGGCTAAGAAGTCAGAGGTAGCCAATGTTTCCTTAAATTCATTTTTAA ACTTACCATTGGTGCATATGTTCCAGATGGCAGATGCTGTCAATAATCTCACCATTGATGACCTTTGTG TATGTAGTTCTTGATCCTATACTGGATAAGCCTGTTTAACTGTATGATGGGTGCTTCCATTGCTT CATAATCTTCATGAAGTGCATGCTTTGCAGCTTTTCACAGTTTATTTGCATTTCTAATGTAGTAATA AAGTAACCAATATAATCATTAA ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTTCGATTCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>Note:</b>	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
<b>RefSeq:</b>	<u><a href="#">NM_021992.3</a></u>
<b>Locus ID:</b>	11013
<b>MW:</b>	16.8