

Product datasheet for **SC205439**

CLUL1 (NM_014410) Human 3' UTR Clone

Product data:

| | |
|--------------------|--|
| Product Type: | 3' UTR Clones |
| Symbol: | CLUL1 |
| Synonyms: | RA337M |
| Mammalian Cell | Neomycin |
| Selection: | |
| Vector: | pMirTarget (PSI00062) |
| ACCN: | NM_014410 |
| Insert Size: | 311 bp |
| Insert Sequence: | <p>>SC205439 3'UTR clone of NM_014410 The sequence shown below is from the reference sequence of NM_014410. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site</p> <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC CATTTTAAGGAACATTTTAAAACTGGGAAGAAGATCTAATGCATCCTATATCCAGTAAGTAGAATTAT CTCTTCATCTGGGACCTGGAAATCCTGAAATAAAAAAGGATAATGCAATAAACACAGTTGCAGGAAAGT ATGTTAGCTATATACTATGAAGTACTCTTAGTTTACTTATGTTGAATGGCTTAGCTATTAATACTCAAA TTGAGTTAAATGAAAATTCCTCCTTAAAAAATCAAACGTAATATGTATTACATTTTCATGGTACATTAG TAGTTCTTTGTATATTGAATAAATACTAAATCACC ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre> |
| Restriction Sites: | SgfI-MluI |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs). |
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |



| | |
|------------------|--|
| Note: | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required. |
| RefSeq: | <u>NM_014410.6</u> |
| Locus ID: | 27098 |
| MW: | 12.5 |