

Product datasheet for **SC205448**

C6orf118 (NM_144980) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: C6orf118 (NM_144980) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: C6orf118
Synonyms: bA85G2.1; dj416F21.2
ACCN: NM_144980
Insert Size: 423 bp
Insert Sequence: >SC205448 3'UTR clone of NM_144980
The sequence shown below is from the reference sequence of NM_144980. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ATTTGTAATAATCAGAGGAAACAGAAGATACGCTCAAGTGGCTGAACAAATGGCCAGTGGAAATTCTTA
ACTTCCATGTATGTAGCATGGATATTTGCAGATGCATATACATTCATAAAGATACATTTTCATTCATCA
GCTTTAATCTCTAGTAAATTAATGTGGTATATGTATCCTACTCAGGATACTGAAATACACACATGCGT
ACACATGATTAGTTTCTGTTGCATACATATTTACAACAATACTAGTTCTAACGTATGTTTTACTTTCA
TTTAATAATGTTCAACTTTGAAAATACTTTTGATTATTTGCTAAGAATAAGGAATTCAGCGCTGCTTG
AATACAGCTCAAATTTGTTTATTTGCCAACACAAGTTTACATATATGCATTTAATAAAATTTTATAGTG
ATTTAATAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_144980.4](#)



[View online »](#)

Locus ID: 168090

MW: 16.4