

Product datasheet for **SC205456**

ZNF226 (NM_001032374) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ZNF226
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_001032374
Insert Size:	426 bp
Insert Sequence:	>SC205456 3'UTR clone of NM_001032374 The sequence shown below is from the reference sequence of NM_001032374. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTTTTGCTCTATGACCTGGCTCAAACCTAACTTGGATTGGAAGTTAGAAGAAATGTTGGAAGTCATT
ATATATGAAGAAATGTTGGAAGGACTCATATATGCATACATTCCTTGAGTGACTATGAATGACTGCCGG
GCAGTAACTTCTGGGCTGTGGTTGTAACCTGTGAGCACTACAAAATGTTTTCTTATTGATACCATAT
TATGGTAGGAAAGACATGGAATAAAAAATTTAGATAGTATGTCAGTAGTTGTGTTTTAAATGGGTTTC
ATTAGTGCTTAGCAATTGGGAGCTTGGTGGACCATCTCTTGGTTTTGGACCATCTCTTGGTTTCTGTCA
GTATGTAAACCAGAACTTCAAATGTGTCACAAAAGATGAGCAGAACTATCCCGAGGTTTATTAAAGTC
TTTTACTCTGTC
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG

```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001032374.2</u>
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	7769
MW:	16.7