

Product datasheet for **SC205457**

ZNF226 (NM_015919) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: ZNF226

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PS100062)

ACCN: NM_015919

Insert Size: 426 bp

Insert Sequence: >SC205457 3'UTR clone of NM_015919
The sequence shown below is from the reference sequence of NM_015919. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTTTTGCTCTATGACCTGGCTCAAACCTAACTTGGATTGGAAGTTAGAAGAAATGTTGGAAGTCATT
ATATATGAAGAAATGTTGGAAGGACTCATATATGCATACATTCCTTGAGTGACTATGAATGACTGCCGG
GCAGTAACTTCTGGGCTGTGGTTGTAACCTGTGAGCACTACAAAATGTTTTCTTATTGATACCATAT
TATGGTAGGAAAGACATGGAATAAAAAATTTAGATAGTATGTCAGTAGTTGTGTTTTAAATGGGTTTC
ATTAGTGCTTAGCAATTGGGAGCTTGGTGGACCATCTCTTGGTTTTGGACCATCTCTGGTTTCTGTCA
GTATGTAAACCAGAACTTCAAATGTGTCACAAAAGATGAGCAGAACTATCCCGAGGTTCAATAAAGTC
TTTTACTCTGTC
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_015919.4</u>
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	7769
MW:	16.7