

Product datasheet for **SC205508**

DTX2 (NM_001102596) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	DTX2 (NM_001102596) Human 3' UTR Clone
Symbol:	DTX2
Synonyms:	RNF58
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001102596
Insert Size:	420 bp
Insert Sequence:	<p>>SC205508 3'UTR clone of NM_001102596 The sequence shown below is from the reference sequence of NM_001102596. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site</p> <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GTGACCGAGGACTGCCTGGAGCAGCAGTGACCTCGCACCCAGCAGCCCCGCTCTGGTGGCCACCCCG CTGCCCATGGCTGGCTGGGTGGCCAGGCAGGAAGTGCCAGCCCGAGAGGCTGGGAGGTTTGTGAGG GTGTGGGGTGTGCCCCACCTGAAGCCGGGGCTCCCCCTGCCTGCCTCTCTCTCCTCCTCCCTCTGGGA ATTGGGCAGCCCTGGGCAGTTGTACTCATGGGGGCTTAGGATGCAGCTACCTCAGTGCGCAGGGCCCGT CTGTCTCTGGGGCTGTCTCGGGCCCGGGTGTCTCGGGGCTGGTGTGGGGCGAGTAGAGACTTCCCC AGCCTGGACGGGCGTGGGTTCTGGGTCAGCTTCTTTTACCTCAATTTTGTGGCAATAAATGCTCTATA GCCAAA ACGCGTAAGCGGCCCGGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online »](#)

RefSeq: [NM_001102596.2](#)

Summary: DTX2 functions as an E3 ubiquitin ligase (Takeyama et al., 2003 [PubMed 12670957]).
[supplied by OMIM, Nov 2009]

Locus ID: 113878

MW: 14.7