

Product datasheet for SC205511

DTX2 (NM_001102594) Human 3' UTR Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | 3' UTR Clones |
| Product Name: | DTX2 (NM_001102594) Human 3' UTR Clone |
| Symbol: | DTX2 |
| Synonyms: | RNF58 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pMirTarget (PS100062) |
| ACCN: | NM_001102594 |
| Insert Size: | 420 bp |
| Insert Sequence: | <p>>SC205511 3'UTR clone of NM_001102594 The sequence shown below is from the reference sequence of NM_001102594. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site</p> <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GTGACCGAGGACTGCCTGGAGCAGCAGTGACCTCGCACCCCAGCACGCCCGCCTCTGGTGGCCACCCCG CTGCCCATGGCTGGCTGGTGGCCAGGCAGGAAGTGCCAGCCCAGAGGCTGGGAGGTTTGTGAGG GTGTGGGGTGTGCCACCTGAAGCCGGGGCTCCCCCTGCCTGCCTCTCTCTCCTCCCTCTGGGA ATTGGGCAGCCCTGGGCAGTTGTACTCATGGGGGCTTAGGATGCAGCTACCTCAGTGCGCAGGGCCCGT CTGTCTCTGGGGCTGCTTCGGGCCCGGGTGCCTCGGGGCTGGTGTGGGGCGAGTAGAGACTTCCCC AGCCTGGACGGGCGTGGGTTCTGGGTCAGCTTCTTTTACCTCAATTTTGTGGCAATAAATGCTCTATA GCCAAA ACGCGTAAGCGGCCCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre> |
| Restriction Sites: | Sgfl-MluI |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs). |
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |



[View online »](#)

RefSeq: [NM_001102594.3](#)

Summary: DTX2 functions as an E3 ubiquitin ligase (Takeyama et al., 2003 [PubMed 12670957]).
[supplied by OMIM, Nov 2009]

Locus ID: 113878

MW: 14.7