

Product datasheet for **SC205518**

SNAPC5 (NM_006049) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: SNAPC5

Synonyms: SNAP19

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_006049

Insert Size: 1021 bp

Insert Sequence: >SC205518 3'UTR clone of NM_006049
The sequence shown below is from the reference sequence of NM_006049. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GAGGAAGAGGAAGAAGAATCAGATTCCAAAGGAGGGGAAAGCCAGGCTGGTGAATGCACATCTGAA
CAGAGACAGGCATGGGAAGGAGGCCAAAAGCGACATGGCAAACCACGGGCAGGGGACTCGCTCTTTG
TTGCTTCCCAACACTTAGACGCCAGCAGCATGGGTTGGTGTGCTGGGCTGGTTAAGGCCGATGGTGA
GCAGAGCCAACCTGCAAAATCCACTTCCTCAGCATCAGATCTCTTGATAAAAGCATGAACCTGCATGTA
AGGAAACGAGAGGACGTGGTGTAAAAAGAAACAAAAACAAACCCACCAGGAAGGAAAATCCAAGA
GTTTCAGGTCTAGAGCTTGTGCTGCACTTGAAAGAGCACCTGGCACCTGCCCTGAGCAACAGTGGGG
CCAGCTGCAGGCCAAAGACCTGCTCCACCAGTCACAGCAGCTGTGAGGCCTTGGGCAAGTCACCTCTT
GGAGCCTACTCGTGCAGTATGGAGACAGTAACCTGCTTAAATGCCTCACAGGCTTGTAGAAGGGAC
AAAGAAAATGATGGAAGTATGATATATATACTATATACACCGTCTATGGCTGTCCAGATATTACA
GTATTACTTAGCTATTTCCACTGCCACCAACGAGACAGCAGGCGGAGAGTGACTCTGGAAGGAGGCTG
AATGTAAGATGAGGCAGTCAGTATTAGGAAAGGCCGCATCCCCATGAAGGGGGTGGGCTGGCCTTTTC
TTTCCCATGGCCCCCTCAGGCATTTCTCTGGCATATGCTTACTCAGGCTGCGGGGGTCACCCCAAGCTT
TCCCTTAAACTGATGGGAGAGCAAGGAGACTGGTGTGTTGTTATTCTGCAGCCCAAGGCTCAGCTGGAC
CACAGTGGGCAGGGAATGAATGGGCTCTGGAAGTACTGCACAGATGACTGACAGAGAAGAACAATGA
TAAACAGGTACTGTAAGAATCCGAATAAATAGACTCACCATGAGTTGCTTCAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_006049.4</u>
Summary:	This gene encodes a subunit of the small nuclear RNA (snRNA)-activating protein complex that plays a role in the transcription of snRNA genes. This complex binds to the promoters of snRNA genes transcribed by either RNA polymerase II or III and recruits other regulatory factors to activate snRNA gene transcription. The encoded protein may play a role in stabilizing this complex. A pseudogene of this gene has been identified on chromosome 6. [provided by RefSeq, Jul 2016]
Locus ID:	10302
MW:	37.8