

Product datasheet for **SC205519**

NTF2 (NUTF2) (NM_005796) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	NTF2 (NUTF2) (NM_005796) Human 3' UTR Clone
Symbol:	NTF2
Synonyms:	NTF-2; NTF2; PP15
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_005796
Insert Size:	1684 bp



[View online »](#)

Insert Sequence: >SC205519 3'UTR clone of NM_005796
The sequence shown below is from the reference sequence of NM_005796. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTCAGGCTCGCCCTGCACAACCTTGGCTGACCTCCTCTCAGCTAGGCACTCACGCTGTTTCCTCCTCCC
TCTCTTCCCAATACTATTCCCACTCCTCCAGATGCTCCAAATATCATGCACAAATGAGCAGGGCCGCG
GTGGGAGTGGGCGCAGTGCCTGCTGCCACTGAGGTGTTGTGCATGATGTTTGGATGCTAGACTAGTTG
CATCTGACGGGAGAAGTTTGTGTGTACCAGCGCATGCCTTGAAAGACTTAAGTAATGAAAAGGTTG
TCCTTTTTTTTTTTTTTTTTTTTTTAATCTACTGACAAGTTGCTCTAGTAACCCAAAGAAGTGAAGGAG
AAAGCAGCTGCCTCACC GCCCAGACATTGATTTGTTTCAGATGTTTCAATGCCTCATGATACAATAAAC
CACAAAAATTTCTTAACAGTTTAAATTGTTTAAATTAGTTTACTAGTTGGCTGGGCATCAGAAGCTAC
CCAGACCGTGTCTCTCTCATGTTTACCCTCCCACTTCTGCCCAACTCTGAATACCTGTTGCAGGG
AGAAACTGCTGAAGCAGCACTCCAGCTAGCCTCTCAGAAGCCGTGTAGGGCAGTTCTCTCCCACTCCA
CAGGGCGACTGGTCTTTGAAAGCAGGCAGGAAGAGCTCAGGTCTGGCCTGAAGCTTCTCTCACTGAAG
CCTGTGCAGTCAGTTCTGATACCTCCTGTGACTTCTGCTCTGGGTAGGTTAGGGCCCTAGTAGCAGGA
GAGATGATCCTGAATCCTGTTGAAGCTGGAGAGGCCTTGCAAAATGGCTCATCACGTTCAAGGCCCTCC
GGGCTGAGTTGTGAGCAGTATCAAGGGAGGGGCTGCTCTATCCCAGAAGGATCAGGATCATATCCAG
GATGCCCCACATACACCAAGCCAGGCAGAGGGCAGCTCAGCTCCTGTCCCATCTGCTTTGGATATCTTT
ACCCAAAGGCAGGTAACCCGAAGAGCCAGCCTCCACTGCCACAGAGCCAGGCCAGTTGTGTTGGAGT
ATAGGTCAGGAGCTGTGGAAGGAGGCAGTCTGTGAGGGACTCATGCTTTAGGAGTCTCACCCCTCAGA
CTGCTGCAGGACATTGCCAGGCCTCTCTCCACTTCTCCTCAGCATACAGACTTCATGCTATCTTCCA
ATTCCGGGAGTCTTAGCTATTAGGTCAGTTTCTGCTTCTCCATTTGGGGACAAAGGCCTTGCCAGT
ACAAATCTAGCCCTTGTCCACAGACTTCTGGATGGTATAAACCTAGTGGCAATGTAGCAACCATAGG
CTAGAACCAAAACCAAGATTTGGGTGAGTGCCTGTTAAGGGTTTAGGATTGGTAAGGACACCACAGC
TAAATCTGACATGTAAAAGGATACCTTCCCCTGTCCCACTACGGGTGGAGGCTAAGGACCTCTCAGA
TCCACAGATTGCCTGGTGACATTGGGCACAGGGCTGGATGGTGTGCTGTTTGGAGGACAGCATGCTC
AGGAGATTTAGTGGGAATGATCTTAGGCTGATGTATTTGGGGAAGTGGGGTGGGGAAGAGGTAG
GTGCTGTGCTGCTCAACCACTTGTGTAAATCTGCAATACGGCTTCTCCATGTACCTGTGCCTGAA
CTGTCTGCAATAAAGAAGATTTGTAAG
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_005796.3](#)

Summary: This gene encodes a cytosolic factor that facilitates protein transport into the nucleus. The encoded protein is required for nuclear import of the small Ras-like GTPase, Ran which is involved in numerous cellular processes. This protein also interacts with the nuclear pore complex glycoprotein p62. [provided by RefSeq, Apr 2016]

Locus ID: 10204

MW: 62.5