

## Product datasheet for **SC205572**

### ERH (NM\_004450) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** ERH (NM\_004450) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** ERH  
**Synonyms:** DROER  
**ACCN:** NM\_004450  
**Insert Size:** 440 bp  
**Insert Sequence:** >SC205572 3'UTR clone of NM\_004450

The sequence shown below is from the reference sequence of NM\_004450. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CGTCGGCAGGCCAACAGGCTGGGAAAATTGTGTTGGAAGCACTGGGGGGTTGGGTGGGCTTGGGA
ACACAGGTGTGTACAGCGTGTGTAGTGAAGTTTTGTATCATAGTAATCCTGTTCCACTTTGTTATA
CTCTAGCCAAGATTGACTGTATTAGATGAAATGTGAGGATCTTGTCAATCGGAAACCCCGTTACCTC
CTCTTTTTCTTTCTTTCTTTTTTTTTTTTTACTTAAACATTTTTATGATGATTTAGATGGAAGTTGT
TCTTCGCACTTAATGTTGGTCCAGTCTTCAACTGTTTCATATCTACTTTATAACATTCACATACTAA
CCCTTCTTCAAGATGGGGTGGGGGTGGAAATGCAGTTTAGCCATGTCCTCAAGATAAAGCTTGGTAA
AAATAAATAAATGTCCTTTAGTTATA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

**Restriction Sites:** Sgfl-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_004450.3](#)



[View online »](#)

**Summary:** May have a role in the cell cycle.[UniProtKB/Swiss-Prot Function]  
**Locus ID:** 2079  
**MW:** 16.7