

Product datasheet for **SC205584**

C19orf48 (NM_199249) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C19orf48 (NM_199249) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	C19orf48
ACCN:	NM_199249
Insert Size:	421 bp
Insert Sequence:	>SC205584 3'UTR clone of NM_199249 The sequence shown below is from the reference sequence of NM_199249. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_199249.3</u>
Locus ID:	84798

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AGAATCCTGGGGAGCTGCTCCTGGTTTGAGCTGCATTCAAGGAAGTGCGGGACATGGTAGGGGAGGCCAA
AAAGCCTTGGGCACTACCCTCCCTGTGGAGCTGTTCCGGTGTCCGTCGAGCTAGCCACACCCTGACACCA
TGTTCAAGGATACCGGAAGAGAAGGGTGTCTGCCCCCAACCTCCCTGTGGGTGTACTGCCAGATGT
CATGAGGGAAGCAGGCCTTGTGAGTGGACACTGACCATGAGTCCCTGGGGGAGTGATCCCCAGGCAT
CGTGTGCCATGTTGCACTTCTGCCAGGCAGCAGGGTGGGTGGGTACCATGGGTGCCACCCCTCCACC
ACATGGGGCCCCAAAGCACTGCAGGCCAAGCAGGGCAACCCACACCCTTGACATAAAAGCATCTTGAA
GCTTTTA
ACGCGTAAGCGGCCCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```



[View online »](#)

MW: 15.3