

## Product datasheet for **SC205740**

### Chk1 (CHEK1) (NM\_001274) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	Chk1 (CHEK1) (NM_001274) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	CHEK1
Synonyms:	CHK1
ACCN:	NM_001274
Insert Size:	1876 bp



[View online »](#)

**Insert Sequence:** >SC205740 3'UTR clone of NM\_001274  
The sequence shown below is from the reference sequence of NM\_001274. The complete sequence of this clone may contain minor differences, such as SNPs.  
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGCCAGAAGATTTGGCTTCTGCCACATGATCGGACCATCGGCTCTGGGAATCCTGGTGAATATAGTG
CTGCTATGTTGACATTATTCTTCTAGAGAAGATTATCCTGTCCTGCAAACTGCAAAATAGTAGTTCCTG
AAGTGTTCACTTCCCTGTTTATCCAAACATCTTCCAATTTATTTTGTGTTGTTGCGCATACAAATAATAC
CTATATCTTAATTGTAAGCAAACTTTGGGAAAGGATGAATAGAATTCATTTGATTATTTCTTCATGT
GTGTTTAGTATCTGAATTTGAACTCATCTGGTGGAAACCAAGTTTCAGGGGACATGAGTTTTCCAGCT
TTTATACACAGTATCTCATTATTTATCAAAACATTTTGTTAATTCAAAAAGTACATATCCATGTTGA
TTAATTCTAAGATGAACCAATAAAGACATAATCTTGACTTTTGGACAGTAGATTTATCAGTCTGT
GAAGCGAAGCCAGCTTCAAAACATATCCCAAGATTTGTACTTATATTTTCAAAAGGCGCTGGCCAGTT
ATATAAACCTGTTTTGAAATTATAATGATTAATTAATAAATTGCAAGTAGGTGTTTTTCCAGTGTAGTTA
GTAAAAACTTGTATTTTACAGTGTTCATAAACTCTAGTGCTTAACTAACTTTACTCTAAAAATTA
GTTGAACATCTTAAATATTTTCTATATTTTCTACTTTCATAGCCATATTTTAACTTTTCACTTACT
GGTGACCAAGCTTTTAGGTGATAAAGAATAAAGAGGGGAAGGGAAGAGTAAGGAAGCTATAAGAAAAAT
AGATCTGATTCTTTGTTCTTTACCTGTTAGACTTACAAAAAGTTTGTGTTTTCTAATAAAATTTGTATC
AACTTTGGGGCATATTAGGTTGAGGCCTGGCTCCTGCCTGTAGTCCCAGCTACTTAGGAGGCTGAGAG
AGGAGGATCGCGTGAACCTGGAAGTTGAGGCTGTAGTGAGCTATGATTGCACCAGTGCCTCCAGCTT
GGATGACAGAGTAAGACCCTACCTCTAATAAAAAATTTTAAAAATGTAACATTATAAAATTAATCAG
TTATTTAATCTGAAGCCAAGAACATGTAGAATGTTATGATTAGAGTTTATCACATATTAATGTACT
GGCAAATTTGTTACTGGAGTATACCCATAGGAGGAATAAATTCAAACCTGTTTTATTTATTTGAACCT
ATTTACGGTATGCTTAAGAATTGAATCAGTATAAATTTCTCAAATATGGGAGAAATTTTGTCTTGAGAA
TTATCTGAGTCATTAATATTTTCAAAAACAGCTCTCACTGACTTGAACCTCTTCTGTAAGCTCTAACC
TTTTACCTGCTTTACATTTCCACTTGAATGTCTAGTAGGCATCTCTTGACCAAAAACAGCTTTTGATTC
CTGTTCTCCAACCTGTTCTCCTAGTTTTCTCCATCTCAGAAATGTTACTTCTCTGCAAAGTCTTT
CCCTGACTTATCTAAAAATAAACCTCCTCTGTTTGTGTTGGGAATTTGTATAGAATGGTGGGAAAAT
TCAAGTTTCATATTTGGATTAGCTCTGACATTTATTTATCTGAACACTGGTAATTCCTCAGTAAAGAC
ACTGATAATAAGTACTTTTAGAGTATTTTAATCTTTAATGCTTAATGTGTAGGAAGAGTATAGTGT
CCTGTTTTGCACAGAAAGCATTCTGTAATAATAAGTTGCCTAATTTTCTGTAATGTTCAATTATAT
TGTGTGGGAAGGTATTTACTCCTATTATTAATAAATAAATGTGTAATAATTTACTACCTGAAAAAAA
AAAAAAAAAAAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** Sgfl-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_001274.5](#)

**Summary:** The protein encoded by this gene belongs to the Ser/Thr protein kinase family. It is required for checkpoint mediated cell cycle arrest in response to DNA damage or the presence of unreplicated DNA. This protein acts to integrate signals from ATM and ATR, two cell cycle proteins involved in DNA damage responses, that also associate with chromatin in meiotic prophase I. Phosphorylation of CDC25A protein phosphatase by this protein is required for cells to delay cell cycle progression in response to double-strand DNA breaks. Several alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Oct 2011]

**Locus ID:** 1111

**MW:** 72.3