

## Product datasheet for SC205755

### Pregnancy specific beta 1 glycoprotein 11 (PSG11) (NM\_002785) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	Pregnancy specific beta 1 glycoprotein 11 (PSG11) (NM_002785) Human 3' UTR Clone
Symbol:	Pregnancy specific beta 1 glycoprotein 11
Synonyms:	PSBG-11; PSBG-13; PSG13; PSG14
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_002785
Insert Size:	457 bp
Insert Sequence:	<p>&gt;SC205755 3'UTR clone of NM_002785</p> <p>The sequence shown below is from the reference sequence of NM_002785. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p>

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCCGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GGAACCTTTGCTTTCAATAATCCAACGTAGCAGCCGTGATGTCATTTTGTATTTTCAGGAAGACTGGCA
GGAGATTTATGAAAAAGACTATGAAAAGGACTCTTGAATACAAGTTCCTGATAACTTCAAGATCATAACC
ACTGGACTAAGAACTTTCAAAATTTTGATGAACAGGCTGATACCTTCATGAAATTCAGACAAAGAAGA
AAAGAACTCCATTTCATTGGACTAAATAACAAAAGGATAATGTTTTCATAATTTTTTATTGGAAATGT
GCTGATTTTTTGAATGTTTTATCCTCCAGATTTATGAATTTTTTCTTCAGCAATTGGTAAAGTATACT
TTTGTAACAAAAATTGAAACATTTGCTTTGCTCTCTGAGTGCCCCAGAAATTGGGAATCTATTCATGA
ATATTCATATGTTTATGGTAATAAAGTTATTTGCACAAGTTTA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

**RefSeq:** NM\_002785.3

**Summary:** The human pregnancy-specific glycoproteins (PSGs) are a group of molecules that are mainly produced by the placental syncytiotrophoblasts during pregnancy. PSGs comprise a subgroup of the carcinoembryonic antigen (CEA) family, which belongs to the immunoglobulin superfamily. For additional general information about the PSG gene family, see PSG1 (MIM 176390).[supplied by OMIM, Oct 2009]

**Locus ID:** 5680

**MW:** 17.7