

## Product datasheet for **SC205767**

### **PWP2H (PWP2) (NM\_005049) Human 3' UTR Clone**

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** PWP2H (PWP2) (NM\_005049) Human 3' UTR Clone  
**Symbol:** PWP2H  
**Synonyms:** EHOC-17; PWP2H; UTP1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pMirTarget (PS100062)  
**ACCN:** NM\_005049  
**Insert Size:** 421 bp  
**Insert Sequence:** >SC205767 3' UTR clone of NM\_005049  
The sequence shown below is from the reference sequence of NM\_005049. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site  
**Blue**=Stop Codon

CAATTGGCAGAGCTCAGAATTCAAGCGATCGC

AGGAGGCAGAGACTCAGAAGAAGAGATGCTGGCCTAGAGCCAGCCGGTTGCAGCGTTGGATTGTGCCGGC  
TAAGACCTGCCAGGGAGATGGGACCCTTGTGCCACCTGGGCCAGCAAAGAGGAGGGGTCCAGAGAACAGC  
TGAATACTGTCACTAGTGGTAGTGACTTGCTTTTCTGTGCACACATGTAGCCCATCAGGACAGCGAGC  
CGACGGGTCACGCCAGGGGCCGGCAGCAGCAGCTGGCACCTGGCCCCAGGAGCGGGCCGTGTGAACGGTGAT  
GAATGTTGAAAATGCGTCTCAGAGAGGTATTCACATGAACTTTGTATGAGACTTATTTATATCTTTAACA  
TAAAGGTTTGATAAAGAAGCTTAGGGATTAAGAAAAAGATGGAGTTTTCTAACGTGAGGATGAAGTCTACA  
C

ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCG

**Restriction Sites:** SgfI-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online >](#)

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_005049.2](#)

**Locus ID:** 5822