

Product datasheet for **SC205787**

APIP (NM_015957) Human 3' UTR Clone

Product data:

| | |
|------------------|-------------------------------------|
| Product Type: | 3' UTR Clones |
| Product Name: | APIP (NM_015957) Human 3' UTR Clone |
| Vector: | pMirTarget (PS100062) |
| Symbol: | APIP |
| Synonyms: | APIP2; CGI-29; CGI29; hAPIP; MMRP19 |
| ACCN: | NM_015957 |
| Insert Size: | 452 bp |
| Insert Sequence: | >SC205787 3'UTR clone of NM_015957 |

The sequence shown below is from the reference sequence of NM_015957. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTCCCAGTTGGAGAAAATGGAATTGTCTAAGCCAAAAGAAAGTCTAATTATACAGAGATAAAGCTAA
ACGTAATTATTATTTAAATGAAAGCTATTTTTTAAATGAATTGAAATTTTTCATGATGCTACTAATTT
GCCACTAAATACTGCAAATGGTCAACCTGAATCTCTTCTGACATTGGATGTTATTTGCTTATATTCTTA
TAATTTTAAATGAGGGCACAGTGAATGAAATTTTATACTCTATGTTTCTGTTTATTTTAAATCCTT
AACAGCAAATATTTGCCTTTAATTTCTTTTTATATACTCTCAGAGAATTCCTCTTAATTTTAAA
GATGCTGGTGATAATAAATTCATTAGAAAATTCCTCATTGTGGAATGAGCATTCTCTGTTTTAATG
TTGGTGCAGAAAATAAATATGAAACATTAAGTCCAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG
```

| | |
|--------------------|--|
| Restriction Sites: | Sgfl-MluI |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| RefSeq: | NM_015957.4 |



[View online »](#)

Summary: APIP is an APAF1 (MIM 602233)-interacting protein that acts as a negative regulator of ischemic/hypoxic injury (Cho et al., 2004 [PubMed 15262985]).[supplied by OMIM, Dec 2008]

Locus ID: 51074

MW: 17.8