

Product datasheet for **SC205871**

C5orf35 (SETD9) (NM_153706) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C5orf35 (SETD9) (NM_153706) Human 3' UTR Clone
Symbol:	C5orf35
Synonyms:	C5orf35
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_153706
Insert Size:	476 bp
Insert Sequence:	>SC205871 3'UTR clone of NM_153706 The sequence shown below is from the reference sequence of NM_153706. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCCGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TTTTCAACTACTACACAATTGTCAGCTAACTCTGTGAATCAGAAATTATTAGGTTTTCTACTCAGCTA
TTAATTCTAAGTGTGTTTTGTTAATCACTTCTCTGAGTTCTACCTGTAAAACAAATATTTTGAGACTTA
ATTGGAATAGGAATTTCTTATGTTTTGTTGATATTATTTATGTTCCAATTTTCATGATAAAAGCTAC
TTGCTTCATTACAATTTCAAATTTGTAATGCAATTCGAAGTTTAATTGGTTTTTCACCTAAATACACAAT
AGCTTATTAAATTAACCTACTCTTTGAACTTATAATTTCAATTTTCTTTGTTATTTTGTAAGCT
CTGTGGTGGTTTTATAAATTGTAGCCAGGCATCAGCAGTTCTTGTAGTACTGCTGATGTGTACAAACTT
CTGTAAAACTTTATTTTACAGAATTGAGTAAAAATACCTATTGTGTTGCCATGAGTAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

RefSeq: NM_153706.4

Locus ID: 133383

MW: 18.6