

Product datasheet for **SC205873**

LDHAL6B (NM_033195) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: LDHAL6B (NM_033195) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: LDHAL6B
Synonyms: LDH6B; LDHAL6; LDHL
ACCN: NM_033195
Insert Size: 532 bp
Insert Sequence: >SC205873 3'UTR clone of NM_033195
The sequence shown below is from the reference sequence of NM_033195. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TGGGAAATTCAGAATAAGCTTAAGCTTTAAAGTTGCCTAAAACCTACCATTCCGAAATTATTGAAGAGAT
CATAGATACAGGATTATATAACGAAATTTGAATAAACTTGAATTCCTAAAAGATGGAAACAGGAAAGT
AGGTAGAGTGATTTTCTATTTATTTAGTCTCCAGCTCTTTTATTGAGCATCCACGTGCTGGACGATA
CTTATTTACAATTCCTAAGTATTTTGGTACCTCTGATGTAGCAGCACTGCCATGTTATATATATGTA
GTTGGCATTGGTTCCCAAAAAGTAGGATGTAGGTATTTATTGTGTTCTAGAAATCCGACTCTTTTCA
TTAGATATATGCTATTTCTTTCATTCTTCTGCTGGTTTATACCTATGTTTATATATGCTGTAAAAAGT
AGTAGCTTCTTACAATGTAATAAATAATGTACATACAAAAAATGCAGTAGTATATACAATCTTTTG
TTTTGCTTCCTTGATAGTTAATAAATCCGTTTGTGAATCAATAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_033195.3](#)



[View online >](#)

Locus ID: 92483

MW: 21.7