

## Product datasheet for **SC205910**

### ZFP38 (ZSCAN21) (NM\_145914) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	ZFP38 (ZSCAN21) (NM_145914) Human 3' UTR Clone
Symbol:	ZFP38
Synonyms:	NY-REN-21; Zipro1; ZNF38
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_145914
Insert Size:	453 bp
Insert Sequence:	>SC205910 3'UTR clone of NM_145914 The sequence shown below is from the reference sequence of NM_145914. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Blue</b> =Stop Codon <b>Red</b> =Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG  
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC  
GTCCACACTGGAGAGGGAGAAGCACCGTACTTTCAAGCGCTCCTGTTGTTGTCGTTGTTTTAACTTT  
AGAATCTGAAAACCAGAAAGAAGTCTTGTCAATGCAGCAGCATCGATTCCGGTGATAGAGTTTGATCA  
CTCAACATCAGGGGATGCCTGAGGAGTGCAGCTCCACAGCAACATGGCAGGCAGGAGTCTCAGAAG  
GTGTCAGGAGTTCCACACTCGCCAGTTCCTGGAGCAGAGTCCCTTCGCCACACTTAGGGTCCCAGTA  
AGCCATGCCAGCATTACCTTTTGCCTAGTTAAACAGACGTGTATCCAGTCTAGTTAAGGAAGAAACATT  
AAGATTGTTTAAATTTTAAACATATATTCAGAATTTTAAATTTGTAAGAATTGAGCCACATTGAACACA  
ATTGAATGAGATTCAGAATAAACTTATAACATCTTGTA  
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA  
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



RefSeq: [NM\\_145914.3](#)

Summary: Strong transcriptional activator.[UniProtKB/Swiss-Prot Function]

Locus ID: 7589

MW: 16.1