

Product datasheet for **SC205917**

CD32 (FCGR2C) (NM_201563) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CD32 (FCGR2C) (NM_201563) Human 3' UTR Clone
Symbol:	CD32
Synonyms:	CD32; CD32C; CDW32; FCG2; FCRIIC; IGFR2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_201563
Insert Size:	1447 bp



[View online »](#)

Insert Sequence: >SC205917 3'UTR clone of NM_201563
The sequence shown below is from the reference sequence of NM_201563. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGCATCGCC
CCCAACGACCATGTCAACAGTAATAACTAAAGAGTAACGTTATGCCATGTGGTCACACTCTCAGCTTGC
TGAGTGGATGACAAAAAGAGGGGAATTGTTAAAGGAAAAATTTAAATGGAGACTGGAAAAATTCCTGAGC
AAACAAAACCACTGGCCCTTAGAAATAGCTTTAACTTTGCTTAACTACAAACACAAGCAAAACTTCA
CGGGGTCATACTACATACAAGCATAAGCAAACTTAACCTGGATGATTTCTGGTAAATGCTTATGTTAG
AAATAAGACAACCCAGCCAATCACAAGCAGCCTACTAACATATAATTAGGTGACTAGGGACTTTCTAA
GAAGATACCTACCCCAAAAAACAATTATGTAATTGAAAACCCATCGATTGCCTTTATTTTGCTTCCAC
ATTTTCCCAATAAATACTTGCCTGTGACATTTTGCCACTGGAACACTAACTTCATGAATTGCGCCTCA
GATTTTTCTTTAACATCTTTTTTTTTGACAGTCTCAATCTGTTACCCAGGCTGGAGTGTAGTGGTGC
TATCTTGGCTCACTGCAACCCGCTCCAGGTTTAAGCGATTCTCATGCCTCAGCTCCAGTAGCTG
GGATTACAGGCATGTGCCGTATAACAGCTAATTTTTGTATTTTTATTTTTTTTTTAGTAGAGACG
GGGTTTCGAATGTTGGCCAGGCCGATCTCGAAGTTCTGGCCTCTAGCGATCTGCCCGCCTCGGCCCTCC
CAAAGTGCTGGGATGACCAGCATCAGCCCCAATGCCAGCCTCTTTAACATCTTCTTTCTATGCCCTC
TCTGTGGATCCCTACTGCTGTTTCTGCCTTCTCCATGCTGAGAACAAAATCACCTATTCAGTCTTAT
GCAGTCGGAAGCTCCAGAAGAACAAAGAGCCCAATTACCAGAACCACATTAAGTCTCATTGTTTGGC
TTGGGATTTGAGAAGAGAATTAGAGAGGTGAGGATCTGGTATTTCTGGTCTAAATCCCCTTGAGGAA
GACGAAGGGATGCTGCAGTTCCAAAAGAGAAGGACTCTCCAGAGTCATCTACCTGAGTCCCGATGCTC
CCTGTCCCTGAAAACACAGACAATATGGTCCCAATGACTGACTGCACCTTCTGTGCCTCAGCCATTCT
TGACATCAAGAATCTTCTGTTCCACATCCACAGCCAATACAATTAATCAAACCACTGTTATGAAAAG
ATGTAGCAACATGAGAAATGCTTATGTTACAGGTTACATGAGAACAATCATGTAAGTCTATATGACTTC
AGAAATGTTAAATAGACTAACCTCTAACACAAATTAAGTGAATGTTTCAAGGTGATGCAATTATT
GATGACCTATTTATTTTCTATAATGATCATATATTACCTTTGTAATAAAACATTATAACCAAAAC
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_201563.5](#)

Summary: This gene encodes one of three members of a family of low-affinity immunoglobulin gamma Fc receptors found on the surface of many immune response cells. The encoded protein is a transmembrane glycoprotein and may be involved in phagocytosis and clearing of immune complexes. An allelic polymorphism in this gene results in both coding and non-coding variants. [provided by RefSeq, Apr 2012]

Locus ID: 9103

MW: 55.2