

Product datasheet for **SC205976**

PRDM11 (NM_020229) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PRDM11 (NM_020229) Human 3' UTR Clone
Symbol:	PRDM11
Synonyms:	PFM8
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_020229
Insert Size:	474 bp
Insert Sequence:	>SC205976 3'UTR clone of NM_020229 The sequence shown below is from the reference sequence of NM_020229. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG  
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC  
AGAGTACGCAGTGGGCTGTGTGGAGGGTAGCCTAAAGCTCTCTGTGGAACACCTTCCGGGAGACCTG  
AGGAGTGAACGTGGAGGCGGCTACCTCCGTGGTGGGAGCCAGGTCTCAGTGTCTCTGGCAGACCC  
ATCGGCAGCTCTGCCAGGTGCTCCATGTGTTGCCCTTGTATCCTCCTTGTCAATAAAGGAAGTCCGCT  
GCAGAAGGGGTGTGTGCTGTGTTCTTGACCCGTTGCCTTTCTCTGGTACTGGTGTCTTACCCCAAAGCC  
CAATTTCTAAACCCAGTCTTTCTGTCCCCAGTCTCAAGCAGGGTGTCCCACTGGAGAGATCTCTTGG  
CTTCCCTAACTTAGTCCAGGAACACAGCCTTGTTCTTCTTCTCCTGAATCTCTGTCTGCCACACATGG  
TCCCAGTTCCTAGCCTGGAGTTCTAGAAGGATGGAGAGTGAGGGGATCCAGGCCATTCA  
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA  
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



RefSeq: [NM_020229.2](#)

Summary: May be involved in transcription regulation.[UniProtKB/Swiss-Prot Function]

Locus ID: 56981

MW: 16.7