

Product datasheet for **SC205987**

SLC44A4 (NM_025257) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: SLC44A4 (NM_025257) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: SLC44A4
Synonyms: C6orf29; CTL4; DFNA72; hTPPT1; NG22; TPPT
ACCN: NM_025257
Insert Size: 465 bp
Insert Sequence: >SC205987 3'UTR clone of NM_025257

The sequence shown below is from the reference sequence of NM_025257. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG  
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC  
CCCCGGACAACAAGAAGAGGAAGAAGTGA CAGCTCCGGCCCTGATCCAGGACTGCACCCACCCCCAC  
CGTCCAGCCATCCAACCTCACTTCGCCTTACAGGTCTCCATTTTGTGGTAAAAAAGGTTTTAGGCCAG  
GCGCCGTGGCTCACGCCTGTAATCCAACACTTTGAGAGGCTGAGGCGGGCGGATCACCTGAGTCAGGAG  
TTCGAGACCAGCCTGGCCAACATGGTGAACCTCCGTCTCTATTAATAAATACAAAAATTAGCCGAGAGT  
GGTGGCATGCACCTGTCTATCCAGCTACTCGGGAGGCTGAGGCAGGAGAATCGCTTGAACCCGGGAGGC  
AGAGGTTGCAGTGAGCCGAGATCGCGCCACTGCACTCCAACCTGGGTGACAGACTCTGTCTCCAAAAACA  
AAACAAACAAACAAAAAGATTTTATTAAGATATTTTGTAACTCAGTAAA  
ACGCGT AAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA  
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_025257.3](#)



[View online »](#)

Summary: The protein encoded by this gene may be a sodium-dependent transmembrane transport protein involved in the uptake of choline by cholinergic neurons. Defects in this gene can cause sialidosis, a lysosomal storage disease. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2010]

Locus ID: 80736

MW: 17.7