

Product datasheet for **SC206003**

BZW2 (NM_014038) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: BZW2 (NM_014038) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: BZW2
Synonyms: HSPC028; MST017; MSTP017
ACCN: NM_014038
Insert Size: 459 bp
Insert Sequence: >SC206003 3'UTR clone of NM_014038
The sequence shown below is from the reference sequence of NM_014038. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GAATCCGAATCGGAAGGTGAGGAAAATTAATGGCTCAACAAGCACAATACCTAGGTTACCACACACCA
CTTTTGTATTGGGAATGCTGAACCATTTGAGAAGAGAACTTGGCTTCTGTTTTCGCAAAGGAAAAAAA
AAATAGGATAGGCTTCCCTTGTGCAGAGGGAGAAATGGTTTTGTTTTGTTTTGTTTTAAATGGAGCC
CTGAGGCATCAGCTATTATACTTGGACTCTACCTCTCACTCACTATATGCTAACTAAAGCCATTCAA
CAAGGAGTCAAGTAGATCTGAAATTAATACTCAACAGACTCCTCTTTTTTAGCTGTATTTTTAGGT
ACTGTGTGGTGACCGCCCACTGGTGTCTATTACAGGCCACTTTGGTAGTTGTGTATCTGCTCATGTAT
GTGATTTGACAAACCAGTTTTTTAAAATAAATGGCTTTTTAAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_014038.3](#)



[View online »](#)

Summary: May be involved in neuronal differentiation.[UniProtKB/Swiss-Prot Function]
Locus ID: 28969
MW: 17.8