

Product datasheet for **SC206020**

CBWD1 (NM_018491) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: CBWD1 (NM_018491) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: CBWD1
Synonyms: COBP
ACCN: NM_018491
Insert Size: 462 bp
Insert Sequence: >SC206020 3' UTR clone of NM_018491
The sequence shown below is from the reference sequence of NM_018491. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site
Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCA**GCGATCGC**

CAGTGGACAACACGTTTCCAAGAAGATCAAGTTTGTACA**TA**CACTAGAGGCATTTCCTTATCAAAAGGAT
TGGATAATAAAAATAAGTTTCTACTGGGTATATTTCAAGCATTATTTATTACTTTAGTTACGAATTCCA
ATATACTTTAAAATGGTATTTGTTTTACAGCATACATAAAATGTAGCAAATCAGTACTGTAAAACATTTA
ACATTCATACAGTTATATATAATATCCTTTTTTTAAAGAATGGTATTTCAAAAAATATCTTTTGAAAT
TGGCTTTGGAGTTTACATATACTGAACATGAAAGTTTATAAATAATGATGATACAACCTTCAACATTGTCA
TTTTTTCTTAGAACTTCAGCTGATTGCAGAGATATAATGATTACATTGTTATTAATTTTTTTAACACAA
GTAAGTGTACCATTTTATGACATGAAATAAAAGGTTATGAC

ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCG

Restriction Sites: Sgfl-Mlul
OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq: [NM_018491.3](#)



[View online »](#)

Locus ID: 55871