

## Product datasheet for **SC206022**

### CBWD1 (NM\_001145356) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	CBWD1 (NM_001145356) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	CBWD1
Synonyms:	COBP
ACCN:	NM_001145356
Insert Size:	476 bp
Insert Sequence:	>SC206022 3'UTR clone of NM_001145356 The sequence shown below is from the reference sequence of NM_001145356. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Blue</b> =Stop Codon <b>Red</b> =Cloning site
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u><a href="#">NM_001145356.1</a></u>

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CGTTTCCAAGAAGATCAAGTTTGTACATAACACTAGAGGCATTTCTTATCAAAAGGATTGGATAATAAA
AATAAGTTTCTACTGGGTATATTTCAAGCATTTATTTACTTTAGTTACGAATCCAATATACTTTA
AAATGGTATTTGTTTACAGCATACATAAAATGTAGCAAATCAGTACTGTAAAACATTTAACATTCATA
CAGTTATATATAATATCCTTTTTTTTTAAAGAATGGTATTTACAAAAATATCTTTTGAAATTGGCTTTG
GAGTTTACATATACTGAACATGAAAGTTTATAATAATGATGATACAACCTTTCAACATTGTCATTTTTTC
TTAGAATTCAGCTGATTGCAGAGATATAATGATTACATTGTTATTAATTTTTTTAACACAAGTAAGT
GTCACCATTTTATGACATGAAATAAAAGTTATGACTGTTATTGAAAAAAAAAAAAAAAAAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG
```



[View online »](#)

Locus ID: 55871

MW: 19.3