

Product datasheet for **SC206022**

CBWD1 (NM_001145356) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CBWD1 (NM_001145356) Human 3' UTR Clone
Symbol:	CBWD1
Synonyms:	COBP
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001145356
Insert Size:	476 bp
Insert Sequence:	<p>>SC206022 3'UTR clone of NM_001145356</p> <p>The sequence shown below is from the reference sequence of NM_001145356. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC CGTTTCCAAGAAGATCAAGTTTGTACAACACTAGAGGCATTTCTTATCAAAAGGATTGGATAATAAA AATAAGTTTCTACTGGGTATATTTCAAGCATTTATTTATTACTTTAGTTACGAATTCCAATATACTTTA AAATGGTATTTGTTTACAGCATACATAAAATGTAGCAAATCAGTACTGTAAAACATTTAACATTCATA CAGTTATATATAATATCCTTTTTTTTAAAGAATGGTATTTACAAAAATATCTTTTGAAATTGGCCTTG GAGTTTACATATACTGAACATGAAAGTTTATAAATGATGATACAACTTTCAACATTGTCAATTTTTTC TTAGAACTTCAGCTGATTGCAGAGATATAATGATTACATTGTTATTAAATTTTTTAAACACAAGTAAGT GTCACCATTTTATGACATGAAATAAAAGTTATGACTGTTATTGAAAAAAAAAAAAAAAAAAAA ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

RefSeq: NM_001145356.1

Locus ID: 55871

MW: 19.3