

Product datasheet for **SC206023**

ARHGAP9 (NM_001080157) Human 3' UTR Clone

Product data:

| | |
|--------------------|--|
| Product Type: | 3' UTR Clones |
| Product Name: | ARHGAP9 (NM_001080157) Human 3' UTR Clone |
| Vector: | pMirTarget (PS100062) |
| Symbol: | ARHGAP9 |
| Synonyms: | 10C; RGL1 |
| ACCN: | NM_001080157 |
| Insert Size: | 445 bp |
| Insert Sequence: | >SC206023 3'UTR clone of NM_001080157 The sequence shown below is from the reference sequence of NM_001080157. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site |
| Restriction Sites: | Sgfl-MluI |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| RefSeq: | <u>NM_001080157.2</u> |

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTGCCCCATTTCCGTGCTGCCCTTGGGTGATAGCACACTCAGATAAGAATCGCATGACACCCCAACAACC
TGGGAATTGTGTTTGGACCAACCCTGTTTCGGCCAGAGCAGGAGACATCTGACCCAGCAGCCCATGCTC
TCTACCCAGGGCAGCTGGTCCAGCTGATGCTACCAACTTCAACAGCCTCTTCCCCTGATGAGGGAAG
GAAGAAGAGAAAACATATTTCCGGTCATCTCTGGTGGTGAAGGGTGGTGTCTGTTTGGAGATATCC
CTTTAAATCTCCAAATGACTGTCTCTATCTTCATGAGTGTGACTTGAGGTGTTGGGATGGGTGAGGGA
GCTTCTCTAAAGAGGAAAGTGAGTGGATTAACCCCTGCTTCTCTTCTTGTTCCTGTTATCATTCTCC
CCAACATAATAATACATAAAGTGGCATGGTA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```



[View online »](#)

Summary: This gene encodes a member of the Rho-GAP family of GTPase activating proteins. The protein has substantial GAP activity towards several Rho-family GTPases in vitro, converting them to an inactive GDP-bound state. It is implicated in regulating adhesion of hematopoietic cells to the extracellular matrix. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Locus ID: 64333

MW: 16.3