

## Product datasheet for **SC206038**

### **CYHR1 (NM\_138496) Human 3' UTR Clone**

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** CYHR1 (NM\_138496) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** CYHR1  
**Synonyms:** CHRP  
**ACCN:** NM\_138496  
**Insert Size:** 478 bp  
**Insert Sequence:** >SC206038 3'UTR clone of NM\_138496  
The sequence shown below is from the reference sequence of NM\_138496. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG  
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC  
CGGCTCTTCTGTCCAGATACAGAAGTAGGGCGGGCCCTCAGGATGTCGAGGAGCCCACGGGCGGCA  
TCCCAGCACCGCTGCCCTGTCCACCTGGCTGGCAGCTGCTTACAGGACTATCTGATCATTAGCAA  
GGAGGAGAAACAAACGAAGCCAACACAGGGCAAGTCTGCATGCGTGCGGACGGGGCCCCCGCTCCGGC  
TACCCCCCGACCCCTGCCTCCCTCCTTCCGAGGGCCGCCAGAGGCTGGGCTGACCCGAAGAGGAGA  
CGGTGCACCAGGCGCCCCGAGGCTCAGAGACGGTGGCAGCAAGGAGGCCGAGAGGCACAGCGACCCCTGC  
CCCAGCCCTTCTGTGCAGTCAGGCGGGGTGCTGCTCCATCCCTGCGGGTTCCGGCGGGGCGCGGGGGC  
CTTGCTGACATCAGACGGGATATCCGAATATCTGATAGCAATTAAGGAGCCTTGTTCGTA  
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA  
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

**Restriction Sites:** Sgfl-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_138496.2](#)



[View online »](#)

Locus ID: 50626

MW: 17