

Product datasheet for **SC206074**

BAP31 (BCAP31) (NM_001139457) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	BAP31 (BCAP31) (NM_001139457) Human 3' UTR Clone
Symbol:	BAP31
Synonyms:	6C6-AG; BAP31; CDM; DDCH; DXS1357E
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001139457
Insert Size:	479 bp
Insert Sequence:	<p>>SC206074 3'UTR clone of NM_001139457</p> <p>The sequence shown below is from the reference sequence of NM_001139457. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC GATGGTCCCATGGACAAGAAGGAAGAGTAAGGCCTCCTTCTCCCTGCCTGCAGCTGGCTTCCACCT GGCAGTGCCTGCTGCTTCTGAGAGCCCGCCTCTCCCTCCAGTACTTGTGTTTGTGCCCTTCTGCTT CCCCATTCCCTTCCACAGCTCATAGCTCGTCATCTCGGCCCTTGTCCACACTCTCCAAGCACATTACA GGGGACCTGATTGCTACACGTTTCAAGATGCGTTTGTGTCATCCTGCTTGGCCTGGCCAGGCCTGGCAC AGCCTTGGCTTCCACGCCTGAGCGTGGAGAGCACGAGTTAGTTGTAGTCCGGCTTGCAGTGGGGCTGAC TTCTGTGTTGGTTTGGCCCTTTTGTGTTTGGCCCTCTGGGTGTTTCTTTGGTCCCGCAGGAGGGTGGG TGGAGCAGGTGGACTGGAGTTTCTCTTGAGGGCAATAAAAGTTGTGATGGTGTGTACGTGGAAAA ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

RefSeq: [NM_001139457.2](#)

Summary: This gene encodes a member of the B-cell receptor associated protein 31 superfamily. The encoded protein is a multi-pass transmembrane protein of the endoplasmic reticulum that is involved in the anterograde transport of membrane proteins from the endoplasmic reticulum to the Golgi and in caspase 8-mediated apoptosis. Microdeletions in this gene are associated with contiguous ABCD1/DXS1375E deletion syndrome (CADD5), a neonatal disorder. Alternative splicing of this gene results in multiple transcript variants. Two related pseudogenes have been identified on chromosome 16. [provided by RefSeq, Jan 2012]

Locus ID: 10134

MW: 16.9