

Product datasheet for **SC206153**

CIDE B (CIDEB) (NM_014430) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: CIDE B (CIDEB) (NM_014430) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: CIDEB
ACCN: NM_014430
Insert Size: 476 bp
Insert Sequence: >SC206153 3'UTR clone of NM_014430
The sequence shown below is from the reference sequence of NM_014430. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CAGCAGAAGGGCCGCCTCATTCTACTAAGGGGCTCTGAGCTTCTGCCCCAGAATCATTCCAACCGA
CCCACTGCAAAGACTATGACAGCATCAAATTTTCAGGACCTGCAGACAGTACAGGCTAGATAACCCACCC
AATTTCCCACTGTCTCTGATCCCTCGTGACAGAACCTTTTCAGCATAACGCCTCACATCCAAGTCT
ATACCCTTACCTGAAGAATGCTGTTCTTTCTAGCCACCTTTCTGGCCTCCCACTTGCCCTGAAAGGCC
AAGATCAAGATGTCCCCAGGCATCTTGATCCCAGCCTGACTGCTGCTACATCTAATCCCCTACCAATG
CCTCTGTCCCTAAACTCCCAGCATACTGATGACAGCCCTCTCTGACTTTACCTTGAGATCTGTCTTC
ATACCCTTCCCCTCAAATAACAAAAACATTTCCAATAAAAAATCAAATATTTACCACTAA
ACGCGTAAGCGGCCCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: Sgfl-Mlul

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_014430.4](#)

Summary: Activates apoptosis.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Locus ID: 27141

MW: 17.9