

Product datasheet for **SC206156**

KBTBD5 (KLHL40) (NM_152393) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	KBTBD5
Synonyms:	KBTBD5; NEM8; SRYP; SYRP
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_152393
Insert Size:	581 bp
Insert Sequence:	>SC206156 3'UTR clone of NM_152393 The sequence shown below is from the reference sequence of NM_152393. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC CTCAATGTGCTGTGCCTGACTAAGATGTGACACAGCTCAGGCAGACTGAACTAAGCACCCCTCCCATCCT GCGACCCTCACTGGCCTGGCCTTGTGGGGCTCCAGAAAAGAGGCTAGGAGAGGCCAGAGTCTACCTGG ATCCAGTTATGGTGCCTCAGGGGCTGCGTCAGCCAAGGAAAGGGAAGTGCTGCTTAGTCTGGACTTTT GGGCAAGGGTGAGAACTAGAGGCTTCTCCAGTGTGCCATATCCCCCTAGGTTGTCTTGATCCATGAA CCAGAACCACAGGGCGGTATCCAGGCCGTGTGCTGGCCCTGCCCCAGCTAGCTGAGTGTGCTGGCAA AGTTCCCCACAGGACTCAGCCTTCTCGTCTGTCCGATGGGAGCATCCCCATCAAGTGCAGGTACAGTG CAGATATGTCTCCTTCTTTAGGAAGAATAAAGTGCCTTCTGAGCAAGCAGCTCAGGGTCTGCATCAGTG CTCCAAGAGTCAGTGAGCAGGTGTGTTGAAGACTTGGGGCTTCAGTGTTCTGCTGTTGGGCTGCCAAGG TCGATGGTCTCTGCTCTCCATCTGGTGTG ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_152393.4</u>
Summary:	This gene encodes a protein containing a BACK domain, a BTB/POZ domain, and 5 Kelch repeats, however, its exact function is not known. The gene and the multi-domain protein structure are conserved across different taxa, including primates, rodents, chicken and zebrafish. [provided by RefSeq, Dec 2012]
Locus ID:	131377
MW:	20.8