

Product datasheet for **SC206179**

GM CSF Receptor alpha (CSF2RA) (NM_172245) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: GM CSF Receptor alpha (CSF2RA) (NM_172245) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: CSF2RA
Synonyms: alphaGMR; CD116; CDw116; CSF2R; CSF2RAX; CSF2RAY; CSF2RX; CSF2RY; GM-CSF-R-alpha; GMCSFR; GMCSFR-alpha; GMR; GMR-alpha; SMDP4
ACCN: NM_172245
Insert Size: 486 bp
Insert Sequence: >SC206179 3'UTR clone of NM_172245

The sequence shown below is from the reference sequence of NM_172245. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GAGGTCTTGACCGTGAAGGAAATTACTGAGACCCAGAGGGTGTAGGAATGGCATGGACATCTCCGCCT
CCGCGACACGGGGAACTGTTTTCTTGATGATGCTGTGAACCTTTATATCATTTTTCTATGTTTTTATTT
AAAAACATGACATTTGGGGCCAGGCGGGTGGCTCACGCCTGTAATCCCAGCACTTTGGGAGGCCAAGG
CAGGCGGATCACCTGAGGTCAGGAGTTCAAGACCAGCCTGCCAACATGGTGAACCCCATCTGGACTA
AAATGCAGAAATTTACCCAGGCACGGCGGCGGACGCCCATCATCCAGCTACTTGGGAGGCTGAGGCA
GGAGAATTGCTTGAACCCGTGAGGCGGAGTTGTAGTGAGCCAAGATCGCACCATTGCACACCAACCTG
CGTGACAGAGCAAGATTGCATCTCAAAACAAACAATAATAATAATAAAAAACCTGATATTTGGCTG
GGC
ACGCGTAAGCGGCCGCGGCATCTAGATTGAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI
OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online »](#)

RefSeq: [NM_172245.4](#)

Summary: The protein encoded by this gene is the alpha subunit of the heterodimeric receptor for colony stimulating factor 2, a cytokine which controls the production, differentiation, and function of granulocytes and macrophages. The encoded protein is a member of the cytokine family of receptors. This gene is found in the pseudoautosomal region (PAR) of the X and Y chromosomes. Multiple transcript variants encoding different isoforms have been found for this gene, with some of the isoforms being membrane-bound and others being soluble. [provided by RefSeq, Jul 2008]

Locus ID: 1438

MW: 18.4