

Product datasheet for SC206202

C15orf40 (NM_001160113) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: C15orf40 (NM_001160113) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: C15orf40
ACCN: NM_001160113
Insert Size: 589 bp
Insert Sequence: >SC206202 3'UTR clone of NM_001160113
 The sequence shown below is from the reference sequence of NM_001160113. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCCGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TCAAGCTATTCTCCTGCCTCAGCTTCCAGTAGCTGGGATTATAGGCGCATGCCACCATGCCTGGCTA
ATTTTTGTATTTTAAATAGAGACGGAGTTTCACCATGTTGGTCAGGCTGGTCTTGAACCTCCTGACCTCA
AATGATCCACCTGCCTTCTCCTCCGAAGCGCTGGGATTAGAGGTATGAGCCACCATGCCCGGAGGTCT
CTCTCTTTTTAACTTCTTTAATTATTTCTTTTTTTTTTTTGGAGACGGAGTCTCGCTCTGTTGTCCAGC
CTGGAGTGCAGTGGTGAATCTTGACGCACTGCTACCTCTCACTCCTGGATTCAAGCGATTCTCCTGCC
TCAGCCTCCTGAGTAGTTGGGATTACAGGCACCCACCACCAGGCCACCTAATTTTTGTATTTGTAGTA
GAGACAGGGTTTACCATGTTGGACAGGCTGGTCTTGAACGCCTGACCTCAGGTGATTCACCCACCCCG
GCCTCCCAAAGTCTGGGATTACAGGCGTGAGCCACTGCGCCTGGCCTTATTTTATTTCTTTTTTGA
TGGAGTCTTGCTCTGTTGCCAGGCTGGAATGCAGT
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 µg dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001160113.2](#)



[View online »](#)

Locus ID: 123207

MW: 21.1