

Product datasheet for **SC206205**

TPRA1 (NM_016372) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	TPRA1
Synonyms:	GPR175; TMEM227; TPRA40
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_016372
Insert Size:	490 bp

Insert Sequence: >SC206205 3'UTR clone of NM_016372
The sequence shown below is from the reference sequence of NM_016372. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGCGAGCGCTGGAAGGCCATCAATGCCGAGGGCAGCTGCCAGGGCCTGTGGAGGACAGGCCAGAGAGG
AGGCCAGCAGGCCAGAGTCCCAGGGGAGGAGGACCAGGTCAAGGGACGTTCTGTGGGCAGTAGCCCT
GTGTGGCCCTGTTCCACCATGAGTCTGGAGGCCACCTCCCTGGGGCTCCCAATCCCTTTGCCATC
TCTGCTCTCACTGGGGACCTCCTCCCTTCCACCTGCTCTCATACTGCTCAGTGACATGGCCAGGC
TTTCCTTCCAGGGCCATGCTTGGCAAGGTTGGCTGAGGGCACCTCCTTCTCTGCACCTTGGCACGAG
GGCAGGGCTGGCTCTCCAATGCCTCCATCCCATCCCATGGTGCTTTGGCCTCCTCAAAGCATCCACC
ATGGTGGATGGACTGAAGTGTATATTTTCTTGATCTATTTTAAATAAAAAGGAAAGGAGCAGAAG
GTGAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_016372.2</u>
Locus ID:	131601
MW:	17.4