

Product datasheet for **SC206222**

BPOZ (ABTB1) (NM_172027) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	BPOZ (ABTB1) (NM_172027) Human 3' UTR Clone
Symbol:	BPOZ
Synonyms:	BPOZ; BTB3; BTBD21; EF1ABP; PP2259
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_172027
Insert Size:	481 bp
Insert Sequence:	>SC206222 3'UTR clone of NM_172027 The sequence shown below is from the reference sequence of NM_172027. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CTGCTCGTGTCCATCGGTCTGGACTGTAGCCCCCTGGCTGGGCAGCCCCAGGGGCCAGGAGCTCTCTT
GGAGACAAGCATGTGTATGCGTTTGTGTGCAGCTCTTCTTCTGCTCCCTGCACATTGAGGGCTTCATG
GGGGGTGCGAGGGGCTCAGTGGGGCTTCTTCCCTCCATGAGCCTGGAGACCCAGGGGAGGATCCAT
TTGGGATGAGCCCCCTCCCCCAATGCACAAGCCAGCCCCAAGACCCTGGGGGTGGACACCACTCAGG
GAAACCTGGGGTGGGGGTGGGCTTTGGTCTTAGCACTTTCTTCTCCAGATCCCCCTACCCACCCAG
TCCCAATCCAGTCTCTGGCCCTTGCCCTAGCCCTGAATTGCTTCTTAAGCTGGTGTTCATGCACA
GGGCCATTGAGGAAGGGCTGGGGGAGTGTGTGGCAATAAAGCTTGAAGGCACCGTGGGAGCATGA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

RefSeq: NM_172027.3

Summary: This gene encodes a protein with an ankyrin repeat region and two BTB/POZ domains, which are thought to be involved in protein-protein interactions. Expression of this gene is activated by the phosphatase and tensin homolog, a tumor suppressor. Alternate splicing results in three transcript variants. [provided by RefSeq, Mar 2010]

Locus ID: 80325

MW: 17.1