

Product datasheet for **SC206224**

CLPTM1 (NM_001294) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: CLPTM1

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PS100062)

ACCN: NM_001294

Insert Size: 474 bp

Insert Sequence: >SC206224 3'UTR clone of NM_001294
The sequence shown below is from the reference sequence of NM_001294. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAGCCAGCAGAGGACAAGAAAAAGGATTAGTCGAGACTGGTCCTCACCTGCTCCGGCTCCTGGCGACCA
CTACCCTGCGTCCCGGCCCTCGCCTCCCTCCCTGTCGCCCTTCCCTGGACAGATCAGGCCGGGG
CGGTGGGAGGCCCGCTCAGGTGAGGCCAGCGTGTGATGTAGGGGCCGGGGCAGGCCAGGTTTGT
TGTGGAGGCGCTGTCTGTCCCTCTGTCCCTCTGTGTTCCAGCCATCTCGCCCTGCCAGCCAGCACCA
CTGGGAATCATGGTGAAGCTGATGCAGCGTTGCCGAGGGGTGGGTTGGCGGGGGTGGGGCCGGGCC
CCCTACGGGATGCCACGGCCGTTTCATCATCTTGTCCCTCGTCCCTACCACACTCCCTCCTAGAC
CGCCGCCCTTTAACACAGTCTGGATTTAATAAAATTCATATGGGTGTTAACTTAACTCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001294.4</u>
Summary:	May play a role in T-cell development.[UniProtKB/Swiss-Prot Function]
Locus ID:	1209
MW:	16.9