

Product datasheet for **SC206244**

ZNF580 (NM_001163423) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ZNF580 (NM_001163423) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ZNF580
ACCN: NM_001163423
Insert Size: 473 bp
Insert Sequence: >SC206244 3'UTR clone of NM_001163423
The sequence shown below is from the reference sequence of NM_001163423. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GAGCTGGCGCAGCACGTGCGCCTCCACTTAGCTCGAGACCCGGCCTGTGCTGCCCTGCCCGTCTCAGGG
CCACCAAGTCTGACCCACACAGCGTCACTCACTCCACACACACCCCTGGCTCTGCTGAGGTTACTGC
CTTACCCTGGGCCTCAGTTTCCCCACCTTCCAAAGGGAGGAGCATCATTCTTCCCTACCCCTTTCTA
GCTGTGTGATGTAGACCAAAGTCGTTGCCCTCCCTGGGCCTGGGAACCGTCCGAACTGGGTTCCAGT
CCAGCTGTGCTGTGAGCCTGTGCAAGTGACATGACCTCTCTAAACCTTGGTTTTCTGCTCTCTGGAG
CGGTGAACCGGTGGTTGTCTGCGGGGAAGAGATGATAAAGAGCACGGGCACGGTCTGGTTCATTTCTGT
ATCTACCCCTTCCGCCACGCCCCCGACCCTTTGCTCAATAAACATTCCGCACTCCA
ACGCGTAAGCGGCCCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001163423.2](#)



[View online »](#)

Summary: Involved in the regulation of endothelial cell proliferation and migration. Mediates H(2)O(2)-induced leukocyte chemotaxis by elevating interleukin-8 production and may play a role in inflammation. May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Locus ID: 51157

MW: 17