

## Product datasheet for **SC206314**

### ZNF799 (NM\_001080821) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF799 (NM_001080821) Human 3' UTR Clone
Symbol:	ZNF799
Synonyms:	HIT-40; ZNF842
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001080821
Insert Size:	480 bp
Insert Sequence:	>SC206314 3'UTR clone of NM_001080821 The sequence shown below is from the reference sequence of NM_001080821. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Blue</b> =Stop Codon <b>Red</b> =Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCCGAAAGATCGCCGTG  
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC  
TCCTTGATAGACATAAAAAGACTCACTAGCATTCTCTCTAAATGTATGGAATGTGGGAAAGCATTAT  
TAATTTTATTTTCATTTCAGATACTGTACGAAACACATTGGAGATAGACCCTGTGAATGTAAGCACTTG  
GTAACACCTTAAGTAATTTCCAGTTTCTTTCCAGTACAGTCATCCCTTGATACCTGCTGGGTATTGGTTC  
CAGCACTCCGTGAGCCATGTCCAGTCCCTTTTATAAAATGACATATTTGTATGTAACCTACCCACATCC  
TCTTGTATACTCTCAATTATGTCTAGATTACTTAAAATACCTCATGCATTGTAAGCTATGCAAAATAG  
TTGTTCTATTGTATTGTTTAGGGACTCATGATAAGGAAAAGAGTCTATGTATTTTCAGTATAGATGCAG  
TAATTGTCAGCCTATCAACATAGTATAGTCAGCCAGAACATTAAGTTTCTTTGTTTCAACTCTCA  
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA  
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



RefSeq: [NM\\_001080821.3](#)

Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Locus ID: 90576

MW: 18.6