

Product datasheet for **SC206327**

TRIM71 (NM_001039111) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TRIM71 (NM_001039111) Human 3' UTR Clone
Symbol:	TRIM71
Synonyms:	HYDCC1; LIN-41; LIN41
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001039111
Insert Size:	2000 bp



[View online »](#)

Summary: The protein encoded by this gene is an E3 ubiquitin-protein ligase that binds with miRNAs and maintains the growth and upkeep of embryonic stem cells. This gene also is involved in the G1-S phase transition of the cell cycle. [provided by RefSeq, Dec 2015]

Locus ID: 131405

MW: 75.6