

Product datasheet for **SC206366**

PRADC1 (NM_032319) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PRADC1 (NM_032319) Human 3' UTR Clone
Symbol:	PRADC1
Synonyms:	C2orf7; PAP21
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_032319
Insert Size:	474 bp
Insert Sequence:	>SC206366 3'UTR clone of NM_032319 The sequence shown below is from the reference sequence of NM_032319. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTGCTGCAACCGCCTGGACCTTCTGGTAGAAGAGTTTGTCCCACATTCCAGCCATAAGTGACTCTGAG
CTGGGAAGGGGAAACCCAGGAATTTTGCTACTTGGAAATTGGAGATAGCATCTGGGGACAAGTGGAGCC
AGGTAGAGGAAAAGGGTTTGGGCGTTGCTAGGCTGAAAGGGAAGCCACACCACTGGCCTTCCTTCCCC
AGGGCCCCAAGGGTGTCTCATGCTACAAGAAGAGGCAAGAGACAGGCCCCAGGGCTTCTGGCTAGAAC
CCGAAACAAAAGGAGCTGAAGGCAGGTGGCCTGAGAGCCATCTGTGACCTGTCACTCACTGGCTCC
AGCCTCCCTACCCAGGGTCTCTGCACAGTGACCTTACAGCAGTTGTTGGAGTGGTTTAAAGAGCTGG
TGTTTGGGGACTCAATAAACCTCACTGACTTTTTAGCAATAAAGCTTCTCATCAGGGTT
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



RefSeq: [NM_032319.3](#)

Summary: Plays a role in the modulation of physical activity and adiposity.[UniProtKB/Swiss-Prot Function]

Locus ID: 84279

MW: 17.2