

Product datasheet for **SC206376**

ZNF227 (NM_182490) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ZNF227 (NM_182490) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ZNF227
ACCN: NM_182490
Insert Size: 468 bp
Insert Sequence: >SC206376 3'UTR clone of NM_182490
The sequence shown below is from the reference sequence of NM_182490. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CAGAAAGTCCATACTGGTAAAAAGCTTTAGAAATGAGAAATGTGTTACCACTTTTGTCTGAATGCACA
TCTTCAAGTTTTGGCTAGTCCATGCTGGTGGTAAACCCTGTAAACTACTGAGAGTGGAAGGGGTTT
GTTCACACTTGGAACTTTCTAACAAATCCATCAAGATGATAACACAGAACCATGAACAGGAATAGAAC
TCGATTTAGGGGAGAAATAGGGCTGGTGGCTCTCTGGTAAGATCTAGTTAATAAATGATCACCTT
TCATTGTGAATATATGCCTGAAGATAATGTGTGGAAGGATATTTGCCATATGCTAACTGGTTTTTGGC
CAGGGAGAGTTTTGGGTATTATCCCTTTCTTAATTTTCATTTTATACTTACAGTGATCATTATTTT
CATAAAAGCTGTAAGCTATGAAAAATGAATAAAATTAATAAAATTTTCTGTA
ACGCGTAAGCGGCCCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_182490.3](#)

Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Locus ID: 7770

MW: 18.4