

Product datasheet for **SC206437**

ZNF804A (NM_194250) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: ZNF804A

Synonyms: C2orf10

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_194250

Insert Size: 496 bp

Insert Sequence: >SC206437 3'UTR clone of NM_194250
The sequence shown below is from the reference sequence of NM_194250. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ACTGTCATCCCTTTGCAACCTCTCTTCTAGTCATCACCATAATGGGAAAAAATACTCTTGAAAACT
ATTGCTATATGCGTTAAGTGTTTCATCTATGTGGGTACATGGCTATTTAACTGGTGGAATAAACTGGCC
GATACATGGCGTCATTGGTTTGAAATCATTACTGTAAGTGCAATGATGCAATAAAATCCCTAAGTTTC
TGATATATAATATTATTAAAGCACTGAATAGTTTGAAAATCAATACAATATATGCTATATATTTAAATG
ATGTCTTAAGAGTATGTATAATGTACATAAAATATATTTATAGTACTCTAATTTATGTTGTAAGTATG
CTCCCTTGGTTTTCTTAATCTTTGTGATTTGCACCTATTTAATGTTTAGACAAAGCTGATGGCACTA
TGTTTTGTATCATGTTCTTGAAACTGTAAATTCAGTAAAAATATCTCTTGCAATAAATTTTGTAA
CTATTTAAGTAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_194250.2</u>
Summary:	The protein encoded by this gene is a zinc finger binding protein. Polymorphisms in this gene, especially rs1344706, are thought to confer increased susceptibility to schizophrenia, bipolar disorder, and heroin addiciton. [provided by RefSeq, Nov 2015]
Locus ID:	91752
MW:	19.3